

**HMG
PLUS**

TECHNICAL SPECIFICATIONS

DISPLAYER	LCD screen
POWER	400Ma Lithum battery
WORKING CONDITON	Temperature 5-40 C(40 104 F), Relative humidity 15 80%RH, Atomsphere Presure: 86-106Kpa
TRANSPORTATION AND STORAGE TEMPERATURE/ HUMIDITY	Temperature 20-55 C(4-131 F), Humidity: 10-93%RH
AUTO SHUT OFF	≤70S
WEIGHT	309g
DIMENSION	37mm x 125mm x 94mm
DIAMETER OF CUFF	22~32cm
ATTACHMENTS	Cuff, 10 test strips, 10 lancets,lithium battery, lancing device, user manuallinclud qe pass, warranty card), giftbox
APPLICATION	Measure human being s blood pressure and blood glucose
LIFETIME	5 Years

MULTIFUNCTIONAL MONITOR

ADVANCED MONITOR TEST BOTH
(BLOOD PRESSURE & BLOOD GLUCOSE)



HOW TO TEST BLOOD GLUCOSE

1. INSERT THE STRIP, MAKE SURE THAT THE CODE APPEARED ON THE METER IS THE SAME AS THE STRIP'S CODE.
2. TAKE THE BLOOD SAMPLE FROM EITHER THE LEFT OR THE RIGHT SIDE OF THE REACTION ZONE.
3. THE TEST RESULT WILL SHOW UP IN 8 SECONDS.
4. THIS METER SET SHOULD BE USED WITH THE MATCHING STRIPS.

HOW TO TEST BLOOD PRESSURE

1. THE LEFT ARM INTO THE CUFF, CUFF AIR PIPE TOWARD THE INSIDE OF THE PALM OF YOUR HAND.
2. WHEN WEARING THE TUBE OUTLET IN THE ARM CUFF EDGE SHOULD BE MIDDLE. FROM THE ELBOW 1-2 CM, AIR TUBE POSITION SHOULD BE ON THE MIDDLE FINGER EXTENDED LINE.
3. WHEN MEASURING, KEEP BODY STRAIGHT, RELAX AND QUIET, KEEP YOUR HEART AT THE SAME HEIGHT WITH THE CUFF CENTER, DON'T MOVE YOUR BODY OR SPEAK, WAIT QUIETLY FOR THE MEASUREMENT TO END.

MEASUREMENT TECHNOLOGY MEASUREMENT RANGE

ACCURACY

COMPRESSING AEROFLUXUS MEMORIES

BLOOD PRESSURE MONITORING

Oscillographic testing
Pressure: 0kPa-37.3kPa (0mmHg-280mmHg)
Pulse value: (30-180) beats/minute
Pressure: within 0.4kPa (3mmHg)
Pulse value: $\pm 5\%$ of the reading
Automatic pressure pump
Automatic rapid exhaust
2*99 memories

BLOOD GLUCOSE MONITORING

Blood glucose monitoring
20-600mg/DL

Within 15 mg/dL when blood glucose <100mg/dl
Within $\pm 15\%$ when blood glucose ≥ 100 mg/dl

180 memories