



Air Mattress With Pump

HMG-MB-00

Mattress Specification: (Bubble Type)

1. Material: PVC 3PLY
2. Thickness: 0.4-0.3mm
3. Size: 200CM*90CM*7CM(inflated)
4. Cell No.: 135pcs

Pump Specification:

1. Size: 21.5X12.0X8.cm
2. Weight: 1.0kgs
3. Voltage: AC 230/220V/110V
4. Frequency: 50HZ/60HZ
5. Power: 10W
6. Pressure range: 20-15kpa
(150-115mmHg)
7. Air output: 5.5-5.0L/minute
8. Cycle time: 10-8 minutes





HANSA MEDICAL GERMANY GmbH

+49(0)610 350 9988 7 info@hansa-medical.de

Robert – bosch – Strabe 32, 63303 Dreieich,
Frankfurt / Germany