

AUTOCLAVE



HMG-A-217 92L

SPECIFICATIONS

Type: Tabletop / Horizontal
Overall Size : 628(W) × 612(H) x 988(L)mm
Chamber Size & Type: 410(Diameter)mm × 699(Depth)mm
92.0 Liter (Cylindrical Type)

Air Removal :Pre-vacuum & Vacuum Dry
Vacuum System: Oil-less Vacuum Pump, 120Liter/min
Loading Tray: 1 EA
Operating Temperature & Pressure : 121° C(1.1Bar) ~ 135° C(2.1Bar)
Display: Graphic LCD Display
Unit Weight: 136 kg
Power Supply: AC220V/60HZ
Power Consumption: 4,600 Watt
Installation Table: 660(Width)×720(Height)×810 (Length)mm
• Optional
Recording :Thermal Dot printer

FEATURES

HMG-A-217, 92L is designed to tabletop and place compatible installation without supply and drainage of water.

- All processes are controlled by micro processor
- Full Automated progressing of all sterilization processes
- All parameters during sterilization processes are being displayed on LCD Monitor
- Multifunctional program of sterilization
- Exclusive sterilization program (134°C PROGRAM & 121°C PROGRAM)
- Function test program (BOWIE-DICK TEST & VACUUM TEST)
- User program (MANUAL PROGRAM & DRY PROGRAM)



HANSA MEDICAL GERMANY GmbH

+49(0)610 350 9988 7 info@hansa-medical.de

Robert – bosch – Straße 32,63303 Dreieich , Frankfurt / Germany
www.hansa-medical.de