



Fetal Doppler

HMG -333



Fetal Doppler



HMG-333

Low Noise & Clear sound

Louder and better sound than handheld doppler. (5W)

Portability

Rechargeable battery and compact design makes portable use.



Trend data save

Store data up to 4 hours and up to 15 patients.

Large display

Larger display than handheld Doppler enables to show more information on the screen.

Central Monitoring Software

Works on BCM software along with other Fetal monitors.



Technical Specification

Category	Fetal Doppler
Model	HMG-333
Frequency (MHz)	2
Detection Range (BPM)	30 - 240
Audio Output (W)	5
Headphone Output	●
Interchangeable probe	●
Waterproof probe	IP X 7
Display type	3.2" Color
Trend display	●
Trend saving	4 Hours
Patient ID	●
Event Mark	●
Quick Guide	●
Multi Language	●
Battery Type	3.7V Li-ion
Battery Run Time (h)	5
Rechargeable Battery	●
BCM-200 Software	●
BCM-350 Central Monitoring Software	●
Weight (g)	1,100
Dimensions (mm) Main Unit (HxWxD)	90 X 250 X 118
Carrying Pouch	○
Optional Accessories	Clamp